

## BRAND TENDON HOLDING FORCEPS

straight; 2x3 teeth

24.0cm (9.6")



**SKU:** 195-241-24

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Seizing Instruments</a>
<b>Aesculap</b>	BF328
<b>KLS Martin</b>	23-676-24
<b>AUTHORS</b>	<a href="#">BRAND</a>
<b>Medicon</b>	78.10.25
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Seizing Instruments</a>
<b>Aesculap</b>	BF328
<b>KLS Martin</b>	23-676-24

## AUTHORS

Medicon

## SPECIFICATIONS

[BRAND](#)

78.10.25

[straight](#), [toothed](#)