

BRAND TENDON HOLDING FORCEPS

straight; 2x3 teeth

15.0cm (6")



SKU: 195-241-15

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Tendon Seizing Instruments |
| Aesculap | BF324 |
| KLS Martin | 23-676-15 |
| AUTHORS | BRAND |
| Medicon | 78.10.15 |
| SPECIFICATIONS | straight , toothed |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Tendon Seizing Instruments |
| Aesculap | BF324 |
| KLS Martin | 23-676-15 |

AUTHORS

Medicon

SPECIFICATIONS

[BRAND](#)

78.10.15

[straight](#), [toothed](#)