

BRONCHIAL SUCTION TUBE

curved; Ø 3.0mm; 25.0cm

SKU: 086-405-30

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| Aesculap | OP362 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| Aesculap | OP362 |