

BROPHY DRESSING FORCEPS

straight; serrated

20.0cm (8"); width of jaw 1.9 +/-0.2

SKU: 020-094-20
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD246
KLS Martin	12-196-20
AUTHORS	BROPHY
Medicon	06.46.20
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD246
KLS Martin	12-196-20

AUTHORS

[BROPHY](#)

Medicon

06.46.20

SPECIFICATIONS

[serrated](#), [straight](#)