

## BROWN NEEDLE HOLDER

straight; diagonally serrated

18.0cm (7.2")



**SKU:** 050-096-18

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>KLS Martin</b>	20-344-18
<b>AUTHORS</b>	<a href="#">BROWN</a>
<b>Medicon</b>	10.15.18
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>KLS Martin</b>	20-344-18
<b>AUTHORS</b>	<a href="#">BROWN</a>
<b>Medicon</b>	10.15.18

## SPECIFICATIONS

[serrated](#), [straight](#)