

BRÜNINGS EXCHANGEABLE TIP

fig. 3; 4.0mm; spoon

for BRÜNINGS interchangeable tubes

SKU: 131-486-03

Categories: [Open Surgery](#), [Goitre](#), [Laryngology](#), [Mouth](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP197
KLS Martin	39-616-03
AUTHORS	BRÜENINGS
Medicon	73.43.14

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP197
KLS Martin	39-616-03
AUTHORS	BRÜENINGS

Medicon

73.43.14