

BRÜNINGS EXCHANGEABLE TIP

fig. 2; 3.0mm; spoon

for BRÜNINGS interchangeable tubes

SKU: 131-486-02

Categories: [Open Surgery](#), [Goitre](#), [Laryngology](#), [Mouth](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP196
KLS Martin	39-616-02
AUTHORS	BRUENINGS

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP196
KLS Martin	39-616-02
AUTHORS	BRUENINGS