

## BRÜNINGS INTERCHANGEABLE TUBE

straight; 32.0cm (12.8")

**SKU:** 131-474-30

**Categories:** [Open Surgery](#), [Goitre](#), [Laryngology](#), [Mouth](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP162
<b>KLS Martin</b>	39-584-30
<b>AUTHORS</b>	<a href="#">BRUENINGS</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP162
<b>KLS Martin</b>	39-584-30
<b>AUTHORS</b>	<a href="#">BRUENINGS</a>