

BRÜNINGS INTERCHANGEABLE TUBE

straight; 41.0cm (16.4")

SKU: 131-474-40

Categories: [Goitre](#), [Laryngology](#), [Mouth](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP164
KLS Martin	39-584-40
AUTHORS	BRUENINGS

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP164
KLS Martin	39-584-40
AUTHORS	BRUENINGS