

BRUNNER INTESTINAL CLAMP

angled; serrated longitudinal

24.0cm (9.6")



SKU: 160-240-24

Categories: [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA191
KLS Martin	30-265-24
AUTHORS	BRUNNER
Medicon	50.29.24
SPECIFICATIONS	angled , serrated longitudinal



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA191
KLS Martin	30-265-24

AUTHORS

[BRUNNER](#)

Medicon

50.29.24

SPECIFICATIONS

[angled](#), [serrated longitudinal](#)