

BRUNNER LIGATURE CONDUCTOR

strong curved


30.0cm (12")

SKU: 056-500-02  
 Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

|                |                                     |
|----------------|-------------------------------------|
| UNIT           | <a href="#">pcs.</a>                |
| TYPE           | <a href="#">Ligature Conductors</a> |
| Aesculap       | BM932                               |
| KLS Martin     | 20-864-30                           |
| AUTHORS        | <a href="#">BRUNNER</a>             |
| Medicon        | 13.18.52                            |
| SPECIFICATIONS | <a href="#">curved</a>              |



ZUSÄTZLICHE INFORMATIONEN

|            |                                     |
|------------|-------------------------------------|
| UNIT       | <a href="#">pcs.</a>                |
| TYPE       | <a href="#">Ligature Conductors</a> |
| Aesculap   | BM932                               |
| KLS Martin | 20-864-30                           |

|                       |                         |
|-----------------------|-------------------------|
| <b>AUTHORS</b>        | <a href="#">BRUNNER</a> |
| <b>Medicon</b>        | 13.18.52                |
| <b>SPECIFICATIONS</b> | <a href="#">curved</a>  |