

BRUNNER LIGATURE CONDUCTOR

right


30.0cm (12")

SKU: 056-500-03
Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM933
AUTHORS	BRUNNER
Medicon	13.18.53



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM933
AUTHORS	BRUNNER
Medicon	13.18.53