

BULLDOG CLAMP ARTERIA; TITAN

strongly curved; 50mm



SKU: 027-045-50

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

MATERIAL [Titanium](#)

KLS Martin 13-145-50

Medicon 15.07.51

SPECIFICATIONS [curved](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

MATERIAL [Titanium](#)

KLS Martin 13-145-50

Medicon 15.07.51

SPECIFICATIONS [curved](#)