

## BULLDOG CLAMP ARTERIA; TITAN

straight; 35mm



**SKU:** 027-042-35

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>KLS Martin</b>	13-142-35
<b>Medicon</b>	15.04.36
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>KLS Martin</b>	13-142-35
<b>Medicon</b>	15.04.36
<b>SPECIFICATIONS</b>	<a href="#">straight</a>