

BULLDOG CLAMP ARTERIA; TITAN

straight; 50mm



SKU: 027-042-50

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-142-50
Medicon	15.04.51
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-142-50
Medicon	15.04.51
SPECIFICATIONS	straight