

BULLDOG CLAMP VENA; TITAN

straight; 35mm



SKU: 027-032-35

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap FB328

MATERIAL [Titanium](#)

KLS Martin 13-132-35

Medicon 15.04.35

SPECIFICATIONS [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap FB328

MATERIAL [Titanium](#)

KLS Martin 13-132-35

Medicon 15.04.35

SPECIFICATIONS

[straight](#)