

BULLDOG CLAMP VENA; TITAN

straight; 40mm



SKU: 027-032-40

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB330
MATERIAL	Titanium
KLS Martin	13-132-40
Medicon	15.04.40
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB330
MATERIAL	Titanium
KLS Martin	13-132-40
Medicon	15.04.40

SPECIFICATIONS

[straight](#)