

BULLDOG CLAMP VENA; TITAN

slightly curved; 50mm



SKU: 027-033-50

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

MATERIAL [Titanium](#)

KLS Martin 13-133-50

Medicon 15.05.50

SPECIFICATIONS [curved](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

MATERIAL [Titanium](#)

KLS Martin 13-133-50

Medicon 15.05.50

SPECIFICATIONS [curved](#)