

BULLDOG CLAMP VENA; TITAN

straight; 30mm



SKU: 027-032-30

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-132-30
Medicon	15.04.30
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
MATERIAL	Titanium
KLS Martin	13-132-30
Medicon	15.04.30
SPECIFICATIONS	straight