

## BUNNELL TENDON STRIPPER

fig. 4; Ø 6.0mm

15.0cm (6")

**SKU:** 195-268-04

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Stripper</a>
<b>Aesculap</b>	FO934
<b>KLS Martin</b>	23-652-04
<b>AUTHORS</b>	<a href="#">BUNNELL</a>
<b>Medicon</b>	78.20.06

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Stripper</a>
<b>Aesculap</b>	FO934
<b>KLS Martin</b>	23-652-04
<b>AUTHORS</b>	<a href="#">BUNNELL</a>

**Medicon**

78.20.06