

BUNNELL TENDON STRIPPERS

Ø 5.0mm; straight

23.0cm (9.2")

SKU: 195-260-23

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
Aesculap	FO940
KLS Martin	23-664-23
AUTHORS	BUNNELL
Medicon	78.26.63
SPECIFICATIONS	straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
Aesculap	FO940
KLS Martin	23-664-23

AUTHORS

[BUNNELL](#)

Medicon

78.26.63

SPECIFICATIONS

[straight](#)