

CASPAR DISC RONGEUR

angled up; 3 x 12mm; toothed

14.0cm (5.6")



SKU: 240-513-03

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR
Medicon	57.23.40
SPECIFICATIONS	angled up , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR
Medicon	57.23.40
SPECIFICATIONS	angled up , toothed