

CASPAR DISC RONGEUR

angled up; 4 x 14mm; toothed

14.0cm (5.6")



SKU: 240-513-04

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR
Medicon	57.23.41
SPECIFICATIONS	angled up , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR
Medicon	57.23.41
SPECIFICATIONS	angled up , toothed