

## CASPAR MICRO CURETTE

3.0mm; serrated; angeled

27.0cm (10 3/4")




**SKU:** 243-801-03

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hypophysectomy Instruments</a>
<b>AUTHORS</b>	<a href="#">CASPAR</a>
<b>Medicon</b>	57.47.43



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hypophysectomy Instruments</a>
<b>AUTHORS</b>	<a href="#">CASPAR</a>
<b>Medicon</b>	57.47.43