

CASTANEDA NEONATAL VASCULAR CLAMP

50° angled; graduated; atraumatic

13.0cm (5.1")



SKU: 142-735-02

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vascular Clamps](#)

KLS Martin 24-593-02

AUTHORS [CASTANEDA](#)

SPECIFICATIONS [angled](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vascular Clamps](#)

KLS Martin 24-593-02

AUTHORS [CASTANEDA](#)

SPECIFICATIONS [angled](#)