

CASTANEDA NEONATAL VASCULAR CLAMP

30° angled; graduated; atraumatic

13.0cm (5.1")



SKU: 142-735-01

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
KLS Martin	24-593-01
AUTHORS	CASTANEDA
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
KLS Martin	24-593-01
AUTHORS	CASTANEDA
SPECIFICATIONS	angled