

## CATHETER FITTING; CONICAL

Ø 2.5mm; 6.5mm



**SKU:** 087-919-01

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stopcocks and Connecting Pieces</a>
<b>KLS Martin</b>	18-825-02

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Stopcocks and Connecting Pieces</a>
<b>KLS Martin</b>	18-825-02