

CHAPUT SEIZING FORCEPS

straight; 2X3 teeth

13.0cm (5.2")



SKU: 040-370-13

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tissue Seizing Forceps
Aesculap	BF300
KLS Martin	13-256-13
AUTHORS	CHAPUT
Medicon	17.30.14
SPECIFICATIONS	straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tissue Seizing Forceps
Aesculap	BF300
KLS Martin	13-256-13

AUTHORS

[CHAPUT](#)

Medicon

17.30.14

SPECIFICATIONS

[straight](#), [toothed](#)