

CHERON DRESSING FORCEPS

shanks angled; serrated

25.0cm (10")


SKU: 030-303-25

Categories: [Open Surgery](#), [Hemostats](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Sponge and Dressing Forceps
Aesculap	BF092
KLS Martin	14-290-25
AUTHORS	CHERON
Medicon	19.28.75
SPECIFICATIONS	angled , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Sponge and Dressing Forceps
Aesculap	BF092
KLS Martin	14-290-25

AUTHORS

[CHERON](#)

Medicon

19.28.75

SPECIFICATIONS

[angled](#), [serrated](#)