

## CLAR-N LED HEAD LAMP 6V

Ø 55mm; with transformer 230V

**SKU:** 100-155-00

**Categories:** [Open Surgery](#), [Diagnostics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Head Lamps</a>
<b>KLS Martin</b>	17-523-04
<b>AUTHORS</b>	<a href="#">CLAR</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Head Lamps</a>
<b>KLS Martin</b>	17-523-04
<b>AUTHORS</b>	<a href="#">CLAR</a>