

CLEANWAVE KERRISON PUNCH; 40° UP CUT

4mm; with ejector; thin footplate

18.0cm (7.2"); 10mm open



SKU: 240-071-04AT-CW

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector