

CLEANWAVE KERRISON PUNCH; 40° UP CUT

5mm; with ejector; thin footplate

18.0cm (7.2"); 10mm open




SKU: 240-071-05AT-CW

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector