

CLEANWAVE KERRISON PUNCH; 40° UP CUT

5mm; with ejector

20.0cm (8"); 10mm open; CC-coated



SKU: 240-072-05ACC-CW

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector