

CLEANWAVE KERRISON PUNCH; 40° UP CUT

2mm, with ejector; thin footplate

23.0cm (9.2"); 10mm open; CC-coated




SKU: 240-073-02ATCC-CW

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 40° up cut, with ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 40° up cut, with ejector |