

CLEANWAVE KERRISON PUNCH; 40° UP CUT

2mm; with ejector; thin footplate

23.0cm (9.2"); 10mm open



SKU: 240-073-02AT-CW

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector