

CLEANWAVE KERRISON PUNCH; 40° UP CUT

5mm, with ejector; thin footplate

23.0cm (9.2"); 10mm open; CC-coated




SKU: 240-073-05ATCC-CW

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , with ejector