

CLIP APPLIER; ROTATABLE; DETACHABLE

for Ethicon clips medium

Ø 5.0mm; 320mm; for Ø 6.5mm trocars



SKU: 612-940-01

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Clip Appliers](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Clip Appliers](#)