

CLIP APPLIER; ROTATABLE; DETACHABLE

for Horizon clips small

Ø 5.0mm; 320mm, for Ø 6.5mm trocars



SKU: 612-946-00

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Clip Appliers](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Clip Appliers](#)