

## CLIP APPLIER; ROTATABLE; LUER LOCK

for SLS clips medium

Ø 10.0mm; 320mm



**SKU:** 612-934-01

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Clip Appliers](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Clip Appliers](#)