

COAKLEY ANTRUM TROCAR

with Luer-Lock tubing connection

Ø 2.3mm; 160mm

SKU: 084-190-16

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Trocars
Aesculap	EJ312
KLS Martin	18-515-16
AUTHORS	COAKLEY
Medicon	09.52.33

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Trocars
Aesculap	EJ312
KLS Martin	18-515-16
AUTHORS	COAKLEY

Medicon

09.52.33