

## COLLIER NEEDLE HOLDER

straight; serrated; fenestrated

12.5cm (5")



**SKU:** 050-022-12

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM207
<b>KLS Martin</b>	20-332-12
<b>AUTHORS</b>	<a href="#">COLLIER</a>
<b>Medicon</b>	10.04.13
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM207
<b>KLS Martin</b>	20-332-12

## AUTHORS

[COLLIER](#)

Medicon

10.04.13

## SPECIFICATIONS

[serrated](#), [straight](#)