


COLLIN DRESSING FORCEPS

straight; serrated

25.0cm (9 3/4")

SKU: 030-300-25
Categories: [Open Surgery](#), [Hemostats](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Sponge and Dressing Forceps
Aesculap	BF091
KLS Martin	14-286-25
AUTHORS	COLLIN
Medicon	19.28.25
SPECIFICATIONS	serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Sponge and Dressing Forceps
Aesculap	BF091
KLS Martin	14-286-25

AUTHORS

[COLLIN](#)

Medicon

19.28.25

SPECIFICATIONS

[serrated](#), [straight](#)