


CONTINUOUS FLOW OUTER SHEATH 28,5FR.

allround perforation; rotatable

for inner sheath 27Fr.; 630-209-28



SKU: 630-209-08

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Sheaths