


## CONTINUOUS FLOW OUTER SHEATH 28,5FR.

allround perforation; rotatable

for inner sheath 27Fr.; 630-209-28



**SKU:** 630-209-08

**Categories:** [Endoscopy](#), [Gynecology](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Sheaths</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Sheaths</a>