

COOLEY ILIAC CLAMP

shanks angled; atraumatic

22.0cm (8.6")



SKU: 142-601-21

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | pcs. |
| TYPE | Vascular Clamps |
| KLS Martin | 24-467-21 |
| AUTHORS | COOLEY |
| SPECIFICATIONS | angled |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---------------------------------|
| UNIT | pcs. |
| TYPE | Vascular Clamps |
| KLS Martin | 24-467-21 |
| AUTHORS | COOLEY |
| SPECIFICATIONS | angled |