

## COOLEY PEDIATRIC VASCULAR CLAMP

fig. 3; 45° angled; atraumatic

14.0cm (5.5")




**SKU:** 142-711-03

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Clamps</a>
<b>KLS Martin</b>	24-573-03
<b>AUTHORS</b>	<a href="#">COOLEY</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Clamps</a>
<b>KLS Martin</b>	24-573-03
<b>AUTHORS</b>	<a href="#">COOLEY</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a>