

COOLEY SUCTION TUBE

curved; with suction control

inner Ø 6.5mm; outer Ø 8.0mm; WL: 20.0cm

SKU: 085-024-02

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| KLS Martin | 18-561-35 |
| AUTHORS | COOLEY |
| Medicon | 09.86.10 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| KLS Martin | 18-561-35 |
| AUTHORS | COOLEY |
| Medicon | 09.86.10 |