

COOLEY SUCTION TUBE

90° curved; suction tip a. shaft removable

outer Ø 6.5mm; AL 13.0cm



SKU: 085-026-90

Categories: [Open Surgery](#), [Puncture / Special Needles](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	COOLEY

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	COOLEY