

## COOLEY SUCTION TUBE

2X curved; suction tip a. shaft removable

outer Ø 8.0mm; WL: 16.0cm

**SKU:** 085-024-01

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-567-33
<b>AUTHORS</b>	<a href="#">COOLEY</a>
<b>Medicon</b>	09.86.20

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-567-33
<b>AUTHORS</b>	<a href="#">COOLEY</a>
<b>Medicon</b>	09.86.20