

COOLEY SUCTION TUBE

2X curved; suction tip a. shaft removable

outer Ø 8.0mm; WL: 16.0cm

SKU: 085-024-01

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
KLS Martin	18-567-33
AUTHORS	COOLEY
Medicon	09.86.20

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
KLS Martin	18-567-33
AUTHORS	COOLEY
Medicon	09.86.20